

PRODUCT	PACK SIZE	% WATER CONTENT	BASE CURVE (mm)	DIAMETER (mm)	POWER RANGE	CYLINDER POWERS	ADD POWERS	AXES	DK (ISO 9913-1 - 1998)	DK/t (@-3.00DS)
<b>AVAIRA</b>										
Avaira®	3	46	8.4 (plus powers), 8.5 (minus powers)	14.2					100	125
<b>BIOFINITY</b>										
Biofinity®	3	48	8.6	14.0					128	160
Biofinity® toric	3	48	8.7	14.5		-0.75 to -2.25DC (0.50D steps)		10° to 180° (10° steps)	128	116
Biofinity® multifocal	3	48	8.6	14.0			+1.00, +1.50, +2.00, +2.50		128	142
<b>BIOMEDICS</b>										
Biomedics® 1 day	30 or 90	52	8.7	14.2					17	24
Biomedics® 1 day extra	30	55	8.6 (minus powers) 8.8 (plus powers)	14.2					19	27
Biomedics® 1 day toric	30	55	8.7	14.5		-0.75, -1.25, -1.75 Call 0870 9000 055 for full range.		20°, 90°, 160°, 180°	19	17
Biomedics® 55 Evolution	3 or 6	55	8.6, 8.9 (minus powers only), 8.8 (plus powers only)	14.2					19	27
Biomedics® toric	3 or 6	55	8.7	14.5		-0.75 to -2.25DC (0.50D steps)		10° to 180° (10° steps)	19	17
<b>PROCLEAR</b>										
Proclear®	3 or 6	62	8.6	14.2					27	42
Proclear® 1 day	30 or 90	60	8.7	14.2					25	28
Proclear® 1 day multifocal	30	60	8.7	14.2			Single power profile up to +2.50		25	28
Proclear® EP	3	60	8.7	14.4			Suitable for up to +1.25 add		25	16
Proclear® multifocal	3	62	8.7	14.4			+1.00 to +2.50 (0.50D steps)		27	17
Proclear® multifocal XR	3	59	8.7	14.4			+1.00 to +4.00 (0.50D steps)		25	Varies by Rx
Proclear® multifocal toric	3	59	8.4, 8.8	14.4		-0.75 to -5.75DC (0.50D steps)	+1.00 to +4.00 (0.50D steps)	5° to 180° (5° steps)	25	Varies by Rx
Proclear® toric	6 pack (8.4mm), 3 or 6 pack (8.8mm)	62	8.4, 8.8	14.4		-0.75 to -2.25DC (0.50D steps)		10° to 180° (10° steps)	27	25
Proclear® toric XR*	3	59	8.8	14.4		-0.75 to -5.75 DC (0.50D steps) See product page.		5° to 180° (5° steps)	25	23

PRODUCT	PACK SIZE	% WATER CONTENT	BASE CURVE (mm)	DIAMETER (mm)	POWER RANGE	CYLINDER POWERS	ADD POWERS	AXES	DK (ISO 9913-1 - 1998)	DK/t (@-3.00DS)	COLOUR
<b>CONVENTIONAL SPHERE</b>											
Lunelle ES 70/Rx	1	70	7.7 to 9.5 in 0.3 steps	13.0 to 15.0 in 0.5 steps	Contact Customer Services to confirm all standard and extended powers				36	18	
Omniflex	1	70	8.4, 8.8	14.3					32	30	
Proclear® sphere	1	59	8.2, 8.5, 8.8	14.2					25	36	
Proclear® sphere Rx	1	59	6.5 to 10.5 in 0.1 steps	12.0 to 17.0 in 0.1 steps					25	36	
Zero 6	1	38	8.4, 8.7, 9.0, 9.3	14.0					10	17	
Zero 6 Rx	1	38	8.1 to 9.6 in 0.3 steps	13.5, 14.0, 14.5, 15.0					10	17	
<b>CONVENTIONAL TORIC</b>											
Lunelle toric	1	70	8.0 to 9.2 (0.3 steps)	14.0, 14.5		-0.50 to -6.00DC (0.25D steps)		1° to 180° (1° steps)	36	18	
Lunelle toric Rx	1	70	8.0 to 9.2 (0.3 steps) 14.3 diameter only. 8.3 to 9.5 (0.3 steps) 15.0 diameter only	14.3, 15.0		-0.50 to -6.00DC (0.20D steps)		1° to 180° (1° steps)	36	24	
Omniflex toric	1	70	8.4, 8.8	14.5		-0.50 to -6.00DC (0.25D steps)		1° to 180° (1° steps)	36	30	
Proclear® Tailor-Made toric	1	59	8.0 to 9.3 (0.1 steps)	13.6 to 15.2 (0.1 steps)		-0.50 to -6.00DC (0.25D steps)		1° to 180° (1° steps)	25	36	
Zero 6 toric	1	38	8.1 to 9.6 (0.3 steps)	14.0		-0.50 to -8.00DC (0.25D steps)		1° to 180° (1° steps)	10	17	
<b>CONVENTIONAL MULTIFOCAL</b>											
Lunelle Variations 70+ UV	1	70	8.0 to 9.2 (0.3 steps)	13.0 to 15.0 (subject to base curve)			Profile A (1.25) Profile B (2.00) Profile C (2.75)		36	18	
<b>COSMETIC</b>											
Expressions® Colors	3	55	8.7	14.4							Dark Blue, Light Blue, Green, Aqua, Grey, Hazel, Brown
Lunelle Colors	1	70	7.4 to 9.5 (0.1 steps)	13.0 to 15.0 (0.1 steps)							Lagoon (blue), Menthe (green), Lemon (yellow), Ambre (brown)
<b>MONTHLY SPHERE</b>											
Frequency® 55	3 or 6	55	8.6, 8.9* (minus powers only), 8.8 (plus powers only)	14.2					19	27	
Frequency® 55 Aspheric	3 or 6	55	8.4* (minus powers only), 8.7	14.4					19	21	
<b>MONTHLY TORIC</b>											
Frequency® XCEL toric	3	55	8.7	14.4		-0.75 to -2.25DC (0.50D steps)		10° to 180° (10° steps)	19	18	
Frequency® XCEL toric XR	3	55	8.7	14.4		-2.75 to -3.75DC (0.50D steps)		5° to 180° (5° steps)	19	18	